

RSNA2016

HITACHI

Inspire the Next

Innovating Healthcare, Embracing the Future



SPEED



SPEED | COMFORT | QUALITY

INSPIRED BY YOUR PATIENTS

Hitachi MR, CT and Ultrasound products are designed to meet the needs of you and your patients with a unique focus on faster exams, greater comfort, and higher image quality. Learn how Hitachi products excel in these attributes.

Hitachi's full product line will be displayed at RSNA in the South Hall, Booth #2511. Please sign onto our website at hitachimed.com/rsna to learn more about our products and to schedule a demo during RSNA.

Magnetic Resonance Imaging	4-5
Computed Tomography	6-7
Ultrasound	8-9
Find us at RSNA	10
RSNA Special Events	11
Give to Charity	12



OASIS

THE HIGHEST FIELD OPEN MRI

Offering unmatched patient comfort, rapid workflow, and the latest clinical capabilities, Oasis 1.2T completes your MRI portfolio ensuring all patient referrals can be handled with the best diagnostic quality available.

SPEED

The open design and the industry's most accommodating patient table help to maximize patient acceptance, streamlining patient positioning for optimal workflow

COMFORT

The unlimited lateral opening provides an unobstructed view ensuring a comfortable experience for even the most challenging patients

QUALITY

Hitachi's unmatched vertical field experience drives the highest field strength open MR to unprecedented image quality



APERTO Lucent PRIME OPEN MRI

APERTO Lucent is an open-sided, mid-field MRI which provides high-field application clarity. With its single pillar design, APERTO Lucent uses open space ergonomics to offer extraordinary comfort to the patient while making examinations completely effortless.

SPEED

Workflow efficiencies at each stage of the exam through his architecture ensure fast exams

COMFORT

The widest patient opening and table of any MRP system offer patient comfort without compromising image quality

QUALITY

Hitachi's continuous improvements along with advanced clinical capabilities provide for high quality, comprehensive imaging



SCENARIA

GREATER SPEED & EASE OF USE

Offering rapid workflow with enhanced comfort and access, Scenaria SE 64/128 is your workhorse CT with greater value to meet your need for lower-dose* with excellent diagnostic quality.

Scenaria's unique All Inclusive pricing bundles equipment and complete support for 5 years. It is an outstanding combination of a premium performance CT and value.

SPEED

0.35 second scan speed with up to 50ips standard reconstruction provides the speed to meet higher workflow demands

COMFORT

Patient centric design combines a wide table-top, larger 75cm aperture, 550lbs bariatric capacity table with 2.0 meter scannable range to ensure enhanced patient access

QUALITY

Hitachi's standard auto lateral-shift table provides 16cm of total shift range to facilitate patient positioning at iso-center for optimum resolution and lower dose

*In clinical use, dose saving features may reduce CT patient dose depending on the clinical task, patient size, anatomical location and clinical practices employed. Consultation with a radiologist and physicist are recommended to determine the appropriate dose needed to obtain diagnostic image quality for a particular clinical task.



Supria

REMARKABLY COMPACT, LOWER DOSE

Supria 16 combines the latest dose saving technologies in a compact design with enhanced patient access to comfortably accommodate a wide range of patients, including bariatric. High quality imaging with automated operator functions speeds workflow.

SPEED

With a streamlined user interface, rapid reconstruction and automation of post-scan tasks, Supria provides accelerated workflow

COMFORT

Supria affords premium patient accommodation features including a 75cm gantry aperture, standard 1.8 meter scannable range table with 500 lb. weight capacity

QUALITY

Dose reduction features included with Supria surpass the XR-29 Standard



ARIETTA

RADIOLOGY CLEARLY DEFINED

The ARIETTA platform offers true shared-service systems. Its wide array of general-purpose and specialty probes allow Radiology departments to support ancillary services – from Radiation Oncology to Neurosurgery – without purchasing multiple specialized systems.

SPEED

Unique dedicated interventional probes along with advanced elastography, contrast ability and image fusion techniques accelerate the path from diagnosis to treatment

COMFORT

Lightweight, hand-contoured probe designs increase scanning comfort

QUALITY

Symphonic Technology carefully preserves and optimizes echo data along the entire image processing chain—from transducer to monitor



3D BREAST ULTRASOUND IN 30 SECONDS

Breast imaging of any kind is often associated with patient anxiety; SOFIA was designed to reduce this through a fast and comfortable exam experience. Powered by Hitachi's premium ARIETTA platform, SOFIA gives providers the peace of mind knowing they are offering their patients the latest in 3D breast ultrasound capabilities.

SPEED

SOFIA captures a 3D volume dataset of an entire breast in only 30 seconds

COMFORT

Its multi-layer memory-foam bed and prone positioning allow patients to remain comfortable and covered during the exam

QUALITY

By presenting radiologists with a single dataset per breast, SOFIA streamlines the interpretation process

HELP HITACHI HELP OTHERS

Just for stopping by the booth and taking a short survey, we'll make a donation to a healthcare charity on your behalf.

CHARITIES OF YOUR CHOICE FOR 2016

Alzheimer's Association

American Cancer Society

RSNA R&E Foundation

St. Jude Children's Research Hospital



COMFORT

HITACHI LOCATION



RSNA EVENTS

SUNDAY, NOVEMBER 27TH, 2016

Opening Session: Digital Revolution in Radiology / 8:30 - 10:15am

Hitachi Exhibit / 10:00am - 5:00pm

MONDAY, NOVEMBER 28TH, 2016

Dr. Pomeranz / MSK Imaging / Hitachi Booth #2511 | 10:30am

Dr. Pomeranz / MSK Imaging / Hitachi Booth #2511 | 1:30pm

Hitachi Exhibit / 10:00am - 5:00pm

TUESDAY, NOVEMBER 29TH, 2016

Dr. Pomeranz / Neuro Imaging / Hitachi Booth #2511 | 10:30am

Dr. Pomeranz / Neuro Imaging / Hitachi Booth #2511 | 1:30pm

Hitachi Exhibit / 10:00am - 5:00pm

WEDNESDAY, NOVEMBER 30TH, 2016

Dr. Pomeranz / Body Imaging / Hitachi Booth #2511 | 10:30am

Hitachi Exhibit / 10:00am - 5:00pm

THURSDAY, DECEMBER 1ST, 2016

Hitachi Exhibit / 10:00am - 2:00pm

STAIRCASE

TRANSPORTATION LOBBY

QUALITY





www.hitachimed.com/RSNA

VISIT US AT RSNA
SOUTH HALL | BOOTH #2511